



Certificate of Analysis

Catalog No.: 1030
Product Name: GEN III MicroPlate™
Lot No.: 3310311
Expiration Date: 31 OCT 2021

Quality Profile

Results


Substrate position check:	Passes test
Substrate response check:	Passes test
ID performance check:	Passes test
Contamination check (No growth at 96 hours):	Passes test

Storage Recommendations

Store between 2 °C and 8 °C

Intended Use

Not for human *in vitro* diagnostic use.


Quality Assurance Approval

11/06/2019
Date

CA1030 Rev B
Effective Date: 11/05/2019



21124 Cabot Blvd
Hayward, CA 94545 USA

Tel: +510 785 2564
Fax: +510 782 4639

www.biolog.com