



# Certificate of Analysis

Catalog No.: 1007  
Product Name: AN MicroPlate™  
Lot No.: 3401231  
Expiration Date: 23 JUL 2021

## Quality Profile

Parameter	Result
Substrate position check:	Passes test
Substrate response check:	Passes test
ID performance check:	Passes test
Contamination check (No growth at 48 hours):	Passes test

Samples of this lot have met performance standards established by Biolog, Inc.

## Storage Recommendations

Store between 2 °C and 8 °C

## Intended Use

Not for human *in vitro* diagnostic use.

Quality Assurance Approval

09 Mar 2020

Date



21124 Cabot Blvd  
Hayward, CA 94545 USA

Tel: +510 785 2564  
Fax: +510 782 4639

[www.biolog.com](http://www.biolog.com)

CA1007  
Rev. C  
Effective date: 12/20/2019