



Certificate of Analysis

Catalog No.: 1030
Product Name: GEN III MicroPlate™
Lot No.: 3609131
Expiration Date: 13 SEP 2024

Quality Profile

Results

Substrate position check:	Passes test
Substrate response check:	Passes test
ID performance check:	Passes test
Contamination check (No growth at 96 hours):	Passes test


Storage Recommendations

Store between 2 °C and 8 °C

Intended Use

Not for human *in vitro* diagnostic use.


Quality Assurance Approval


Date

CA1030 Rev C
Effective Date: 12/20/2019



21124 Cabot Blvd
Hayward, CA 94545 USA

Tel. +510 785 2564
Fax. +510 782 4639

www.biolog.com