



21124 Cabot Blvd.
Hayward, California 94545
510-785-2564 Telephone
800-284-4949 Orders
510-782-4639 Fax
www.biolog.com

CERTIFICATE OF ANALYSIS

Product No.: 3419
Name: GP-ROD SB Turbidity Standard
Lot No.: B092617
Expiration Date: 26 JUN 2019

QUALITY PROFILE

	Results
Standard check: 28.0 ± 0.5 % Transmittance	Passes test


Note: The % transmittance values are obtained from specific reference turbidimeters. These values may differ from turbidimeter to turbidimeter.

STORAGE RECOMMENDATIONS

Store in the dark at 20 °C – 26 °C.

INTENDED USE

NOT FOR HUMAN IN VITRO DIAGNOSTIC USE.



APPROVAL: QUALITY ASSURANCE



DATE