



## SAFETY DATA SHEET

### 1 Product and Company Identification

Product name: **Biolog Redox Dye Mix E (100x)**  
Product number: 74225  
Manufacturer/Supplier: Biolog Inc.  
21124 Cabot Blvd.  
Hayward, CA 94545  
Telephone: 1-510-785-2564  
Information department: Technical services.  
Emergency information: In case of a chemical emergency, spill, fire, exposure, or accident contact Biolog Inc. at 1-510-785-2564 or 1-800-284-4949

### 2 Hazards Identification

Hazard description: This product contains no hazardous constituents, or the concentrations of all chemical constituents are below the regulatory threshold limits described by Occupational Safety Health Administration Hazard Communications Standard 29 CFR 1910.1200 and the European Directive 91/155/EEC, and 93/112/EC.  
Classification system: This classification was made based on information gathered from company and literature data.

NFPA ratings (scale 0-4)



Health = 0  
Fire = 0  
Reactivity = 0

### 3 Composition/Information

Chemical characterization:  
Description: Mixture of the substances listed below with nonhazardous additions.  
Dangerous components: Void

### 4 First aid measures

General information: No special measures required.  
After inhalation: Seek medical treatment in case of complaints.  
After skin contact: Immediately wash with water and soap and rinse thoroughly.  
After eye contact: Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.  
After swallowing: If symptoms persist consult doctor.  
Information for doctor: Show this SDS.

### 5 Fire fighting measures

Suitable extinguishing agents: CO<sub>2</sub>, ABC multipurpose dry chemical or water spray. Fight larger fires with water spray or alcohol resistant foam.  
Protective equipment: No special measures required.

### 6 Accidental release measures

Person-related safety precautions: None required.  
Measures for environmental protection: Wipe up with damp sponge or mop.  
Measures for cleaning/collecting: No special measures required.  
Additional information: No dangerous substances are releases.

### 7 Handling and storage

Handling	
Information on safe handling	No special measures required.
Information about protection against explosions and fires	No special measures required.
Storage	
Requirements to be met by storerooms and receptacles	2- 8 °C in the dark
Information about storage in one common storage facility	No special measures required.

### 8 Exposure controls and personal protection

Components with limit values that require monitoring at the workplace	Contains no substances with occupational exposure limit values.
Personal protective equipment	
General protective and hygienic measures	Use usual precautionary measures for handling chemicals.
Breathing equipment	None required.
Protection of hands	Chemical resistant gloves (i.e. nitrile, or equivalent)
Eye protection	Safety glasses
Body protection	Protective work clothing (lab coat)

### 9 Physical and chemical properties

General Information	
Form	Liquid
Color	Pale Yellow
Odor	Not applicable
Specific gravity	Not determined
Density	Not determined
Evaporation rate	Not determined
Boiling point/ Freezing point	Not determined
Solubility in/ Miscibility with Water	Soluble

### 10 Stability and reactivity

Stability	Stable
Thermal decomposition/conditions to be avoided	No decomposition if used according to specifications.
Hazardous decompositions products	Carbon monoxide, Carbon dioxide, Nitrogen oxides, Hydrogen chloride gas.
Hazardous polymerization	Will not occur.

### 11 Toxicological information

Acute toxicity	
Primary irritant effect	
On the skin	May cause irritation.
On the eye	May cause irritation.
Sensitization	No sensitizing effects known
Additional toxicological information	

### 12 Ecological information

Toxicity	No data available
Persistence and degradability	No data available
Bioaccumulative potential	No data available
Mobility in soil	No data available
Results of PBT and vPvB assessment	No data available
Other adverse effects	No data available

### 13 Disposal considerations

Product	
Recommendation	Smaller quantities can be disposed of with solid waste. Dispose of material in accordance with federal (40 CFR 261.3), state and local requirements. This product is not considered RCRA hazardous waste.
Uncleaned package	
Recommendation	Disposal must be made according to state and federal regulations.
Recommended cleaning agent	Water, if necessary with cleaning agents.

**14 Transportation information**

DOT Considered to be non-hazardous for transport.  
IATA Considered to be non-hazardous for air transport.

**15 Regulatory information**

United States regulatory information  
SARA Listed No

Canada regulatory information

WHMIS classification Classified in accordance with the hazard criteria of the CPR, and the SDS contains all information required by the CPR.

DSL No  
NDSL No

California Proposition 65 – Chemicals known to cause cancer  
None of the ingredients listed.

California Proposition 65 - Chemicals known to cause reproductive toxicity  
None of the ingredients listed.

**16 Other information**

*The information, data, and recommendations contained herein are believed to be accurate. Biolog, Inc. makes no warranty of any kind whatever with respect thereto and disclaims all liability from reliance thereon. We reserve the right to revise this SDS periodically as new information becomes available.*

MS74225 Dye Mix E  
Date: 31 March 2021