

Customer Name: <Last Name>, First Name>
Customer Company: <Customer Name>
Customer Address: <Street>, <City>, <State> <ZIP> <Country>

Run Date: MM/DD/YY HH:MM AM/PM
Sample Number: MYMMDD-NNNN: Customer Sample Code
Library Version: (Bruker) 6.00, (MIDI Labs) 08/18/20

Closest Match: <Genus species>
Confidence Level: Choose an item.
Note: Choose an item.

Top 10 Matches:

Rank	Score	Organism Name	Reference Source
1	Score	Name	Source
2	Score	Name	Source
3	Score	Name	Source
4	Score	Name	Source
5	Score	Name	Source
6	Score	Name	Source
7	Score	Name	Source
8	Score	Name	Source
9	Score	Name	Source
10	Score	Name	Source

Score Value Key	
Range	Confidence Level
2.000 – 3.000	Species
2.000 – 3.000 multiple species	Species, Closely Related
1.700 – 1.999	Genus
0.000 – 1.699	No Match

Created: 10/2/2023 2:26 PM
Reviewer: Joshua West 10/2/2023 2:27 PM