

# CERTIFICATE OF ANALYSIS

**Catalog No.:** 3441  
**Product Name:** 85% Turbidity Standard  
**Lot No.:** M100725  
**Expiration Date:** 07 JUL 2027

---

## Quality Profile

---

### Results

Standard check: 85.0 ± 0.5% transmittance      Passes test

Note: The % transmittance values are obtained from specific reference turbidimeters. These values may differ from turbidimeter to turbidimeter. Reference turbidimeters are calibrated with NIST traceable standards.

---

## Storage Recommendations:

Store in the dark between 20 °C and 26 °C.

---

## Intended Use

Not for human *in vitro* diagnostic use

For Research Use Only

---

  
\_\_\_\_\_  
**Quality Assurance**

10 OCT 2025  
\_\_\_\_\_  
**Date**

**biolog**

21124 Cabot Blvd.  
Hayward, CA 94545

biolog.com  
+1 800 284 4949

CA3441 RevB  
Effective Date: 10/20/2023