

CERTIFICATE OF ANALYSIS

Catalog No.: 74225
Product Name: Redox Dye Mix E (100X)
Lot No.: D1104251
Expiration Date: 04 NOV 2027

Quality Profile

<u>Challenge organisms</u>	<u>Source</u>	<u>Growth Result</u>
----------------------------	---------------	----------------------

No Growth Promotion Required

Final pH: 7.1 ± 0.1 at 25 °C
Color: Light Orange, Clear Liquid
Unit: Bottle
Volume: 20 mL

Samples of this lot have met pH, appearance, and performance standards established by Biolog, Inc. Representative samples were incubated at 26, 30 and 35 °C for 96 hours for contamination testing and found to be acceptable. Results were obtained at the time of release. Certification does not exempt the user from his/her responsibility to assure satisfactory performance at the time of use.

Storage Recommendations

Store between 2 °C and 8 °C

Intended Use

Not for human *in vitro* diagnostic use
For Research Use Only

Quality Assurance

12 NOV 2025

Date

biolog

21124 Cabot Blvd.
Hayward, CA 94545

biolog.com
+1 800 284 4949

CA74225 RevD
Effective Date: 10/20/2023