



21124 Cabot Blvd.
 Hayward, California 94545
 510-785-2564 Telephone
 800-284-4949 Orders
 510-782-4639 Fax
 www.biolog.com

CERTIFICATE OF PERFORMANCE

Product No.: 80102
 Name: RAINBOW® AGAR 0157
 Lot No.: 605191
 Expiration Date: 06/ 2019

QUALITY PROFILE

Test 1. Efficiency of plating compared to TSA+B medium

Rainbow O157 comparable at:
 10 uL of approx. 10³ cells/mL Passes test

Test 2. Colony coloration of 16 bacterial strains

Species tested

E. coli	K. pneumoniae	
E. coli O157 (2 strains)	E. amnigenus	
E. coli O111:H8 (2 strains)	S. sonnei	
E. coli O26:H11 (2 strains)	P. aeruginosa	
E. coli O48:H21	E. faecalis	
H. alvei (2 strains)	S. epidermidis	Passes test

Test 3. Verification of Indole Spot Test

Species tested

E. coli	positive	
E. coli O157:H7	positive	
H. alvei A	negative	
P. aeruginosa	negative	Passes test

Samples of this lot have met performance standards established by Biolog, Inc.

Technician: SR Completion date: 05/11/2016

INTENDED USE

NOT FOR HUMAN IN VITRO DIAGNOSTIC USE.


 APPROVAL: QUALITY ASSURANCE


 DATE