



21124 Cabot Blvd.
Hayward, California 94545
510-785-2564 Telephone
800-284-4949 Orders
510-782-4639 Fax
www.biolog.com

CERTIFICATE OF PERFORMANCE

Product No.: 80102
Name: RAINBOW® AGAR 0157
Lot No.: 712181
Expiration Date: JAN 2021

QUALITY PROFILE

Test 1. Efficiency of plating compared to TSA+B medium

Rainbow O157 comparable at:
10 uL of approx. 10^3 cells/mL Passes test

Test 2. Colony coloration of 16 bacterial strains

Species tested

E. coli	K. pneumoniae
E. coli O157 (2 strains)	E. amnigenus
E. coli O111:H8 (2 strains)	S. sonnei
E. coli O26:H11 (2 strains)	P. aeruginosa
E. coli O48:H21	E. faecalis
H. alvei (2 strains)	S. epidermidis

Passes test

Test 3. Verification of Indole Spot Test

Species tested

E. coli	positive
E. coli O157:H7	positive
H. alvei A	negative
P. aeruginosa	negative

Passes test

Samples of this lot have met performance standards established by Biolog, Inc.

Technician: Sarah Roppiyakuda Completion date: 14 DEC 2017

INTENDED USE

NOT FOR HUMAN IN VITRO DIAGNOSTIC USE.


APPROVAL: QUALITY ASSURANCE

05 JAN 2018
DATE