

Catalog No.: 72501
Product Name: YT Inoculating Fluid
Lot No.: 1810301
Expiration Date: 30 OCT 2020

Quality Profile

Parameter	Result
Appearance:	Passed
Contamination check:	Passed

Samples of this lot have met appearance and sterility standards established by Biolog, Inc. Representative samples were incubated at 26, 30 and 35 °C for 96 hours for contamination testing and found to be acceptable. Results were obtained at time of release. Certification does not exempt the user from his/her responsibility to assure satisfactory performance at the time of use.

Color: Clear, colorless liquid
Unit: Tube
Volume: 12.5 mL

Intended Use

Not for human *in vitro* diagnostic use.



Quality Assurance Approval

07 November 2018
Date