



21124 Cabot Blvd.  
Hayward, CA 94545  
510-785-2564 Telephone  
800-284-4949 Orders  
510-782-4639 Fax  
www.biolog.com

### CERTIFICATE OF PERFORMANCE

**Catalog No.:** 70007  
**Product Name:** BUA™ Dehydrated Media  
**Lot No.:** 905072  
**Expiration Date:** 06 / 2022

#### QUALITY PROFILE

Quality Control Organism	Source (ATCC)	Medium	Incubation	Growth Response
<i>Bifidobacterium breve</i>	15700(der.)	5% Blood	Anaerobic 35 °C	Growth
<i>Clostridium sordellii</i>	9714(der.)	5% Blood	Anaerobic 35 °C	Growth
<i>Lactobacillus casei</i>	393(der.)	5% Blood	Anaerobic 35 °C	Growth
<i>Micromonas micros</i>	33270(der.)	5% Blood	Anaerobic 35 °C	Growth

**Final pH:** 7.2 ± 0.1 at 25 °C, BUA with Blood

Samples of this lot have met pH, appearance, and performance standards established by Biolog, Inc.

**Technician:** Shah Rappato **Completion Date:** 21 JUN 2019

#### INTENDED USE

Not for human *in vitro* diagnostic use.

Shah Rappato  
Quality Assurance Approval

07/03/2019  
Date