



## Certificate of Analysis

Catalog No.: 72301  
Product Name: Biolog IF-M1  
Lot No.: 0610211  
Expiration Date: 06 OCT 2024

---

### Quality Profile

---

#### Results

Color: Clear with a faint yellow tint  
Performance check: Passes test  
Contamination check: Passes test

---

### Storage Recommendations

---

Store between 2 °C to 8 °C

---

### Intended Use

---

Not for human *in vitro* diagnostic use.

  
Quality Assurance Approval

01 Nov 2021  
Date

CA 72301 Rev A  
Effective Date: 01/15/2021



21124 Cabot Blvd  
Hayward, CA 94545 USA

Tel: +510 785 2564  
Fax: +510 782 4639

[www.biolog.com](http://www.biolog.com)