



Certificate of Analysis

Catalog No.: 72301
Product Name: Biolog IF-M1
Lot No.: 3105221
Expiration Date: 31 MAY 2025

Quality Profile

Results

Color: Clear with a faint yellow tint
Performance check: Passes test
Contamination check: Passes test

Storage Recommendations

Store between 2 °C to 8 °C

Intended Use

Not for human *in vitro* diagnostic use.


Quality Assurance Approval

28 June 2022
Date

CA 72301 Rev A
Effective Date: 01/15/2021



21124 Cabot Blvd
Hayward, CA 94545 USA

Tel: +510 785 2564
Fax: +510 782 4639

www.biolog.com