



21124 Cabot Blvd.  
Hayward, CA 94545  
510-785-2564 Telephone  
800-284-4949 Orders  
510-782-4639 Fax  
www.biolog.com

### CERTIFICATE OF PERFORMANCE

**Catalog No.:** 70005  
**Product Name:** BUY™ Dehydrated Media  
**Lot No.:** 2302BH  
**Expiration Date:** 18 MAR 2026

#### QUALITY PROFILE

Quality Control Organism	Source	Incubation	Growth Response
<i>Candida albicans</i>	NCYC 1363	Aerobic 26 °C	Growth
<i>Papiliotrema laurentii</i>	CBS 139	Aerobic 26 °C	Growth
<i>Schwanniomyces polymorphus</i> <i>var polymorphus</i>	NCYC 948	Aerobic 26 °C	Growth
<i>Kluyveromyces marxianus</i>	GSU/AHEARN C90006070	Aerobic 26 °C	Growth
<i>Priceomyces haplophilus</i>	NCYC 549	Aerobic 26 °C	Growth
<i>Leucosporidiella muscorum</i>	NCYC 693	Aerobic 26 °C	Growth

**Final pH:** 5.6 ± 0.6 at 25 °C

Samples of this lot have met pH, appearance, and performance standards established by Biolog, Inc.

**Technician:** Antonio Munguia **Completion Date:** 05/05/2023

#### INTENDED USE

Not for human *in vitro* diagnostic use.

David Pina  
Quality Assurance Approval

05/05/2023  
Date