



21124 Cabot Blvd.
Hayward, CA 94545
510-785-2564 Telephone
800-284-4949 Orders
510-782-4639 Fax
www.biolog.com

CERTIFICATE OF PERFORMANCE

Catalog No.: 70007
Product Name: BUA™ Dehydrated Media
Lot No.: 610191
Expiration Date: 11/2019

QUALITY PROFILE

Quality Control Organism	Source (ATCC)	Medium	Incubation	Growth Response
<i>Bifidobacterium breve</i>	15700(der.)	5% Blood	Anaerobic 35 °C	Growth
<i>Clostridium sordellii</i>	9714(der.)	5% Blood	Anaerobic 35 °C	Growth
<i>Lactobacillus casei</i>	393(der.)	5% Blood	Anaerobic 35 °C	Growth
<i>Micromonas micros</i>	33270(der.)	5% Blood	Anaerobic 35 °C	Growth

Final pH: 7.2 ± 0.1 at 25 °C, BUA with Blood

Samples of this lot have met pH, appearance, and performance standards established by Biolog, Inc.

Technician: PT _____ **Completion Date:** 17OCT2016 _____

INTENDED USE

Not for human *in vitro* diagnostic use.


Quality Assurance Approval

07 NOV 2016
Date