



Certificate of Performance

Catalog No.: 72264
Product Name: PM IF-10a GP Base (1.2X)
Lot No.: 2305251
Expiration Date: 25 MAY 2024

Quality Profile

| <u>Challenge organisms</u> | <u>Source</u> | <u>Growth Result</u> |
|------------------------------|---------------|----------------------|
| NO GROWTH PROMOTION REQUIRED | | |

Final pH: 7.1 ± 0.2 at 25 °C.
Color: Light yellow, clear liquid
Unit: Bottle
Volume: 125 mL

Samples of this lot have met pH, appearance, and performance standards established by Biolog, Inc. Representative samples were incubated at 26, 30 and 35 °C for 96 hours for contamination testing and found to be acceptable. Results were obtained at time of release. Certification does not exempt the user from his/her responsibility to assure satisfactory performance at the time of use.

Intended Use

For research use only.


Quality Assurance Approval

06/13/2023
Date



21124 Cabot Blvd
Hayward, CA 94545 USA

Tel: +510 785 2564
Fax: +510 782 4639

www.biolog.com

CA72264
Rev C
Effective Date: 12/20/2019