



Certificate of Performance

Catalog No.: 72501
Product Name: YT Inoculating Fluid
Lot No.: 2307211
Expiration Date: 21 JUL 2025

Quality Profile

| Parameter | Result |
|----------------------|-------------|
| Appearance: | Passes test |
| Contamination check: | Passes test |

Samples of this lot have met appearance and sterility standards established by Biolog, Inc. Representative samples were incubated at 26, 30 and 35 °C for 96 hours for contamination testing and found to be acceptable. Results were obtained at time of release. Certification does not exempt the user from his/her responsibility to assure satisfactory performance at the time of use.

Color: Clear, colorless liquid
Unit: Tube
Volume: 12.5 mL

Intended Use

Not for human *in vitro* diagnostic use.


Quality Assurance Approval

07/26/2023
Date



21124 Cabot Blvd
Hayward, CA 94545 USA

Tel: +510 785 2564
Fax: +510 782 4639

www.biolog.com

CA72501
Rev. B

Effective Date: 12/20/2019