

CERTIFICATE OF ANALYSIS

Catalog No.: 72301
Product Name: Biolog IF-M1
Lot No.: 2009241
Expiration Date: 20 SEP 2027

Quality Profile

Parameter	Results
Color:	Clear with a faint yellow tint
Performance check:	Passes test
Contamination check	Passes test

Storage Recommendations

Store between 2 °C to 8 °C

Intended Use

Not for human *in vitro* diagnostic use

For Research Use Only


Quality Assurance

04 OCT 2024
Date

biolog

21124 Cabot Blvd.
Hayward, CA 94545

biolog.com
+1 800 284 4949

CA72301 RevB
Effective Date: 10/20/2023