

# CERTIFICATE OF ANALYSIS

**Catalog No.:** 72303  
**Product Name:** Biolog Mitochondrial Assay Solution (2x)  
**Lot No.:** 0403251  
**Expiration Date:** 03 APR 2027

---

## Quality Profile

---

No growth promotion required.

**Color:** Colorless, clear liquid  
**Unit:** Bottle  
**Volume:** 100 mL

---

## Storage Recommendations

Store between 2 °C to 8 °C

---

## Intended Use

Not for human *in vitro* diagnostic use

For Research Use Only

---

  
\_\_\_\_\_  
Quality Assurance

18 APR 2025  
Date